

## LABORATORY SUPPLY REQUISITION

DIRECT ORDERS TO GML CLIENT SERVICES PHONE: 800-695-6491

FAX: 570-271-5660

DATE:	# ORDERED	#ISSUED	LAB USE
ORDERED BY NAME/SITE ADDRESS AND PHONE:			ONLY
ALIQUOT TUBES / CAPS (100/PKG)			0150
BAG, SPEC TRANSPORT (6"x9") (100/PKG) RM TEMP REFRIG FROZEN			0150
CONTAINER, STERILE SPECIMEN (90 ML.) (EACH OR 100/PKG)			0090
CULTURE eSWAB - (BACTERIAL, FUNGAL & GC) (EACH)			0050
CULTURE SWAB - GREEN (BORDETELLA CULTURE) (EACH)			0170
CULTURE SWAB - BLUE (GRP A PROBE & RAPID STREP) (EACH/PKG)			0170
CULTURE SWAB - RED (STAPH AUREUS NASAL SCREEN, VRE RECTAL)			0170
CYTOLOGY FIXATIVE or CYTOLOGY PRESERVATIVE TUBES (EACH)			0170
FORMALIN (EACH) OR (32/PKG) 20mL or 40 mL			0090
GEN PROBE GC/CHLAMYDIA COLLECTION KITS			0150
UNISEX SWAB (URETHRAL/ENDOCERVICAL) (EACH/BOX)			0200
URINE COLLECTION KIT (EACH/BOX)			0170
MULT-TEST COLLECTION SWAB (EACH)			0170
LABELS: SPECIMEN IDENTIFICATION			0170
OTHER (SPECIFY):			0150
ECLIPSE NEEDLES (with Barrels): (21 G X 1¼") (22 G X 1¼") (100/Box or Each)			0150
PIPETTE, TRANSFER 7 ML. (EACH OR 500/BOX)			0150
REQUISITIONS:			0150
CLINICAL CUSTOMIZED (EACH) SUPPLY (EACH)			0130
CYTOPATHOLOGY (EACH)  SURGICAL PATHOLOGY (EACH)			0150
SLIDE MAILERS - (20/PKG)			0150
STOOL TRANSPORT VIALS (CULTURE/GIARDIA EIA) (EACH OR 50/BOX)			0150
STOOL TRANSPORT VIALS (O&P) (EACH OR 50/BOX)			0170
SURE PATH ® PAP COLLECTION VIALS & BROOMS (25/PKG)			0170
TUBES, BLOOD COLLECTION (EACH OR 100/PKG)			0150
BLUE, LIGHT (COAGULATION) 4.5 ML.			0130
BLUE, DARK (STERILE) 7 ML. or DK BLUE EDTA 7 ML			0150
GRAY, 3 ML.			0150
GREEN (LITHIUM HEPARIN) 3 ML. or (SODIUM HEPARIN) 3 ML or 10 ML.			0150
LAVENDER (EDTA) 3 ML.			0150
PINK (EDTA - BLOOD BANK) 7 ML.			0150
RED (PLAIN) 7 ML. or GOLD (SERUM SEPARATOR) 4 ML.			0150
TUBE, TRACE ELEMENT FOR LEAD TESTING ONLY, 6 ML			0150
OTHER - SPECIFY:			0150
TUBES, BLOOD COLLECTION - MICROCOLLECTION (EACH OR 20/PKG)			0150
AMBER (BILIRUBINS) RED (PLAIN) GREEN (LITHIUM HEPARIN)			
LAVENDER (EDTA) YELLOW (SERUM SEPARATOR)			0150
URINE ALIQUOT TUBE (24 HOUR), 50 ML. (EACH OR 25/PKG)			0150
URINE CONTAINERS (24 HOUR) (EACH):	1	1	0050
PLAIN ACID ADDED OTHER:			
URINALYSIS TRANSPORT TUBE, 8 ML. (EACH)			0050
URINE CULTURE COLLECTION KIT (CUP AND WIPE) (EACH)			0050
URINE CULTURE TRANSPORT KIT (STRAW AND TUBE TOGETHER) (EA OR BOX OF 50)			0170
URINE TRANSFER STRAW (EACH)			0170
UTM w/FLOCKED SWAB (EACH)			0170
OTHER:			
RECEIVED BY: DATE:	•	•	•